

Baculoviral IAP repeat-containing protein 5

Catalog No: #AP81364



Package Size: #AP81364-1 50ug #AP81364-2 100ug #AP81364-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Baculoviral IAP repeat-containing protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Apoptosis inhibitor survivin
Accession No.	Uniprot ID:Q9JHY7
Uniprot	Q9JHY7
GeneID	64041;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Birc5

Note: This product is for in vitro research use only