

Voltage-dependent P/Q-type calcium channel subunit alpha-1A

Catalog No: #AP81366

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81366-1 50ug #AP81366-2 100ug #AP81366-3 1mg

Description

Product Name	Voltage-dependent P/Q-type calcium channel subunit alpha-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CHCACHA1A, Voltage-gated calcium channel subunit alpha Cav2.1
Accession No.	Uniprot ID:O73705
Uniprot	O73705
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CACNA1A

Note: This product is for in vitro research use only