

## Bactericidal permeability-increasing protein

Catalog No: #AP81368



Package Size: #AP81368-1 50ug #AP81368-2 100ug #AP81368-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bactericidal permeability-increasing protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BPI
Accession No.	Uniprot ID:Q67E05
Uniprot	Q67E05
GeneID	329547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Bpi

Note: This product is for in vitro research use only