

ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1

Catalog No: #AP81371



Package Size: #AP81371-1 50ug #AP81371-2 100ug #AP81371-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2'-phospho-ADP-ribosyl cyclase,2'-phospho-ADP-ribosyl cyclase/2'-phospho-cyclic-ADP-ribose transferase,2'-phospho-cyclic-ADP-ribose transferase,ADP-ribosyl cyclase 1,ADPRC 1,Cyclic ADP-ribose hydrolase 1,cADPr hydrolase 1,I-19,NIM-R5 antigen,CD38
Accession No.	Uniprot ID:P56528
Uniprot	P56528
GeneID	12494;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd38

Note: This product is for in vitro research use only