

CD44 antigen

Catalog No: #AP81373



Package Size: #AP81373-1 50ug #AP81373-2 100ug #AP81373-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD44 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Extracellular matrix receptor-III,ECMR-III,GP90 lymphocyte homing/adhesion receptor,HUTCH-I,Hermes antigen,Hyaluronate receptor,Phagocytic glycoprotein 1,PGP-1,Phagocytic glycoprotein I,PGP-I,CD44
Accession No.	Uniprot ID:Q28284
Uniprot	Q28284
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD44

Note: This product is for in vitro research use only