

Cocaine- and amphetamine-regulated transcript protein

Catalog No: #AP81375

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81375-1 50ug #AP81375-2 100ug #AP81375-3 1mg

Description

Product Name	Cocaine- and amphetamine-regulated transcript protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CART
Accession No.	Uniprot ID:Q16568
Uniprot	Q16568
GeneID	9607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CARTPT

Note: This product is for in vitro research use only