

Calsequestrin-1

Catalog No: #AP81377



Package Size: #AP81377-1 50ug #AP81377-2 100ug #AP81377-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Calsequestrin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Calmitine,Calsequestrin, skeletal muscle isoform,CASQ
Accession No.	Uniprot ID:P31415
Uniprot	P31415
GeneID	844;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CASQ1

Note: This product is for in vitro research use only