

Cell adhesion molecule 3

Catalog No: #AP81385



Package Size: #AP81385-1 50ug #AP81385-2 100ug #AP81385-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cell adhesion molecule 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Brain immunoglobulin receptor, Immunoglobulin superfamily member 4B, IgSF4B, Nectin-like protein 1, NECL-1, Synaptic cell adhesion molecule 3, SynCAM3, TSLC1-like protein 1, TSL1, IGSF4B, NECL1, SYNCAM3, TSLL1, UNQ225/PRO258
Accession No.	Uniprot ID:Q8N126
Uniprot	Q8N126
GeneID	57863;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CADM3

Note: This product is for in vitro research use only