

Calbindin

Catalog No: #AP81388



Package Size: #AP81388-1 50ug #AP81388-2 100ug #AP81388-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Calbindin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Calbindin D28,D-28K,Vitamin D-dependent calcium-binding protein, avian-type,CAB27
Accession No.	Uniprot ID:P05937
Uniprot	P05937
GeneID	793;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CALB1

Note: This product is for in vitro research use only