

ERK1 Mouse Monoclonal Antibody(5E9)

Catalog No: #44201

Package Size: #44201-1 50ul #44201-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ERK1 Mouse Monoclonal Antibody(5E9)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5.00E+09
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	ERK1 Mouse Monoclonal Antibody(5E9) detects endogenous levels of MAPK3
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	ERK1
Other Names	MAPK3
Accession No.	Swiss-Prot:P27361 Gene ID:5595
Uniprot	P27361
GeneID	5595;
Calculated MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using ERK1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ERK1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only