

## Protein chibby homolog 1

Catalog No: #AP81392



Package Size: #AP81392-1 50ug #AP81392-2 100ug #AP81392-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein chibby homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ARPP-binding protein,Cytosolic leucine-rich protein,PIGEA-14,PKD2 interactor, Golgi and endoplasmic reticulum-associated 1,ARB1,C22orf2,CBY,PGEA1,HRIHFB2025
Accession No.	Uniprot ID:Q9Y3M2
Uniprot	Q9Y3M2
GeneID	25776;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CBY1

Note: This product is for in vitro research use only