

Ceramide synthase 1

Catalog No: #AP81394

Package Size: #AP81394-1 50ug #AP81394-2 100ug #AP81394-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ceramide synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CerS1,LAG1 longevity assurance homolog 1,Longevity assurance gene 1 protein homolog 1,Protein UOG-1,LAG1,LASS1,UOG1
Accession No.	Uniprot ID:P27544
Uniprot	P27544
GeneID	10715;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CERS1

Note: This product is for in vitro research use only