

Cadherin-1

Catalog No: #AP81398

Package Size: #AP81398-1 50ug #AP81398-2 100ug #AP81398-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cadherin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ARC-1, Epithelial cadherin, E-cadherin, Uvomorulin, CD324
Accession No.	Uniprot ID:P09803
Uniprot	P09803
GeneID	12550;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdh1

Note: This product is for in vitro research use only