

Carboxypeptidase M

Catalog No: #AP81404



Package Size: #AP81404-1 50ug #AP81404-2 100ug #AP81404-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carboxypeptidase M
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CPM
Accession No.	Uniprot ID:P14384
Uniprot	P14384
GeneID	1368;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CPM

Note: This product is for in vitro research use only