

C-C motif chemokine 15

Catalog No: #AP81407



Package Size: #AP81407-1 50ug #AP81407-2 100ug #AP81407-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chemokine CC-2,HCC-2,Leukotactin-1,LKN-1,MIP-1 delta,Macrophage inflammatory protein 5,MIP-5,Mrp-2b,NCC-3,Small-inducible cytokine A15,MIP5,NCC3,SCYA15
Accession No.	Uniprot ID:Q16663
Uniprot	Q16663
GeneID	6359;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL15

Note: This product is for in vitro research use only