

## C-C motif chemokine 4

Catalog No: #AP81408

Package Size: #AP81408-1 50ug #AP81408-2 100ug #AP81408-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C motif chemokine 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-26 T-lymphocyte-secreted protein,HC21,Lymphocyte activation gene 1 protein,LAG-1,MIP-1-beta(1-69),Macrophage inflammatory protein 1-beta,MIP-1-beta,PAT 744,Protein H400,SIS-gamma,Small-inducible cytokine A4,T-cell activation protein 2,ACT-2,LAG1,MIP1B,SCYA4
Accession No.	Uniprot ID:P13236
Uniprot	P13236
GeneID	388372;6351;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CCL4

Note: This product is for in vitro research use only