

C-C chemokine receptor type 2

Catalog No: #AP81409

Package Size: #AP81409-1 50ug #AP81409-2 100ug #AP81409-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-C CKR-2,JE/FIC receptor,MCP-1 receptor,CD192,Cmkbr2
Accession No.	Uniprot ID:P51683
Uniprot	P51683
GeneID	12772;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccr2

Note: This product is for in vitro research use only