

# Myeloid cell surface antigen CD33

Catalog No: #AP81410



Package Size: #AP81410-1 50ug #AP81410-2 100ug #AP81410-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Myeloid cell surface antigen CD33                                   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Hu  |
| Other Names        | Sialic acid-binding Ig-like lectin 3, Siglec-3, gp67, CD33, SIGLEC3 |
| Accession No.      | Uniprot ID:P20138   |
| Uniprot            | P20138  |
| GeneID             | 945;  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)              |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: CD33

Note: This product is for in vitro research use only