

Carbonic anhydrase 12

Catalog No: #AP81411



Package Size: #AP81411-1 50ug #AP81411-2 100ug #AP81411-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbonic anhydrase 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Carbonate dehydratase XII, Carbonic anhydrase XII, CA-XII
Accession No.	Uniprot ID:O43570
Uniprot	O43570
GeneID	771;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CA12

Note: This product is for in vitro research use only