

Breast cancer type 1 susceptibility protein homolog

Catalog No: #AP81415

Package Size: #AP81415-1 50ug #AP81415-2 100ug #AP81415-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Breast cancer type 1 susceptibility protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	RING-type E3 ubiquitin transferase BRCA1
Accession No.	Uniprot ID:Q864U1
Uniprot	Q864U1
GeneID	353120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: BRCA1

Note: This product is for in vitro research use only