

C-C motif chemokine 28

Catalog No: #AP81418



Package Size: #AP81418-1 50ug #AP81418-2 100ug #AP81418-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 28
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Small-inducible cytokine A28,SCYA28
Accession No.	Uniprot ID:Q68A91
Uniprot	Q68A91
GeneID	448795;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL28

Note: This product is for in vitro research use only