

C-C motif chemokine 5

Catalog No: #AP81419



Package Size: #AP81419-1 50ug #AP81419-2 100ug #AP81419-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MuRantes,SIS-delta,Small-inducible cytokine A5,T-cell-specific protein RANTES,Scya5
Accession No.	Uniprot ID:P30882
Uniprot	P30882
GeneID	20304;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl5

Note: This product is for in vitro research use only