

CD166 antigen

Catalog No: #AP81420



Package Size: #AP81420-1 50ug #AP81420-2 100ug #AP81420-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD166 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Activated leukocyte cell adhesion molecule, CD166
Accession No.	Uniprot ID:Q9BH13
Uniprot	Q9BH13
GeneID	281614;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ALCAM

Note: This product is for in vitro research use only