

Cystathionine beta-synthase

Catalog No: #AP81423



Package Size: #AP81423-1 50ug #AP81423-2 100ug #AP81423-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystathionine beta-synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-thionase,Hemoprotein H-450,Serine sulfhydrase
Accession No.	Uniprot ID:P32232
Uniprot	P32232
GeneID	24250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cbs

Note: This product is for in vitro research use only