

Cyclin-dependent kinase inhibitor 2A, isoform 1

Catalog No: #AP81424



Package Size: #AP81424-1 50ug #AP81424-2 100ug #AP81424-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase inhibitor 2A, isoform 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cyclin-dependent kinase 4 inhibitor A,CDK4I,p16-INK4a,p16,p16-INK4,P16ink4a
Accession No.	Uniprot ID:Q9R0Z3
Uniprot	Q9R0Z3
GeneID	25163;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdkn2a

Note: This product is for in vitro research use only