

CCAAT/enhancer-binding protein beta

Catalog No: #AP81426



Package Size: #AP81426-1 50ug #AP81426-2 100ug #AP81426-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CCAAT/enhancer-binding protein beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C/EBP beta,AGP/EBP,Interleukin-6-dependent-binding protein,IL-6DBP,Liver-enriched transcriptional activator,LAP
Accession No.	Uniprot ID:P28033
Uniprot	P28033
GeneID	12608;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cebpb

Note: This product is for in vitro research use only