

Centromere protein R

Catalog No: #AP81427



Package Size: #AP81427-1 50ug #AP81427-2 100ug #AP81427-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Centromere protein R
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CENP-R,Nuclear receptor-interacting factor 3,Cenpr,Nrif3
Accession No.	Uniprot ID:Q5U1Z7
Uniprot	Q5U1Z7
GeneID	362548;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ltgb3bp

Note: This product is for in vitro research use only