

Protein atonal homolog 1

Catalog No: #AP81428

Package Size: #AP81428-1 50ug #AP81428-2 100ug #AP81428-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Protein atonal homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class A basic helix-loop-helix protein 14,bHLHa14,Helix-loop-helix protein hATH-1,hATH1,ATH1,BHLHA14
Accession No.	Uniprot ID:Q92858
Uniprot	Q92858
GeneID	474;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATOH1

Note: This product is for in vitro research use only