

Ceramide synthase 1

Catalog No: #AP81434



Package Size: #AP81434-1 50ug #AP81434-2 100ug #AP81434-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Ceramide synthase 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | CerS1,LAG1 longevity assurance homolog 1,Longevity assurance gene 1 protein homolog 1,Protein UOG-1,Lass1,Uog-1,Uog1 |
| Accession No. | Uniprot ID:P27545 |
| Uniprot | P27545 |
| GeneID | 93898; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cers1

Note: This product is for in vitro research use only