

B-lymphocyte antigen CD19

Catalog No: #AP81436



Package Size: #AP81436-1 50ug #AP81436-2 100ug #AP81436-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	B-lymphocyte antigen CD19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Differentiation antigen CD19,CD19
Accession No.	Uniprot ID:P25918
Uniprot	P25918
GeneID	12478;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd19

Note: This product is for in vitro research use only