

Cyclin-dependent-like kinase 5

Catalog No: #AP81437



Package Size: #AP81437-1 50ug #AP81437-2 100ug #AP81437-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent-like kinase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cell division protein kinase 5,Proline-directed protein kinase 33 kDa subunit,PDPK,Serine/threonine-protein kinase PSSALRE,Tau protein kinase II catalytic subunit,TPKII catalytic subunit,CDKN5
Accession No.	Uniprot ID:Q02399
Uniprot	Q02399
GeneID	281066;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDK5

Note: This product is for in vitro research use only