

Calpain-1 catalytic subunit

Catalog No: #AP81449



Package Size: #AP81449-1 50ug #AP81449-2 100ug #AP81449-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calpain-1 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Calcium-activated neutral proteinase 1,CANP 1,Calpain mu-type,Calpain-1 large subunit,Micromolar-calpain,muCANP
Accession No.	Uniprot ID:P35750
Uniprot	P35750
GeneID	397027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CAPN1

Note: This product is for in vitro research use only