

## Cathepsin L1

Catalog No: #AP81450



Package Size: #AP81450-1 50ug #AP81450-2 100ug #AP81450-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cathepsin L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cathepsin L,Cyclic protein 2,CP-2,Major excreted protein,MEP,Ctsl1
Accession No.	Uniprot ID:P07154
Uniprot	P07154
GeneID	25697;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ctsl

Note: This product is for in vitro research use only