

T-cell surface antigen CD2

Catalog No: #AP81453



Package Size: #AP81453-1 50ug #AP81453-2 100ug #AP81453-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-cell surface antigen CD2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Erythrocyte receptor,LFA-2,LFA-3 receptor,Rosette receptor,T-cell surface antigen T11/Leu-5,CD2,SRBC
Accession No.	Uniprot ID:P06729
Uniprot	P06729
GeneID	914;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD2

Note: This product is for in vitro research use only