

CD2-associated protein

Catalog No: #AP81454



Package Size: #AP81454-1 50ug #AP81454-2 100ug #AP81454-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD2-associated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mesenchyme-to-epithelium transition protein with SH3 domains 1,METS-1,Mets1
Accession No.	Uniprot ID:Q9JLQ0
Uniprot	Q9JLQ0
GeneID	12488;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd2ap

Note: This product is for in vitro research use only