

## C-C motif chemokine 13

Catalog No: #AP81455



Package Size: #AP81455-1 50ug #AP81455-2 100ug #AP81455-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C motif chemokine 13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CK-beta-10, Monocyte chemoattractant protein 4, Monocyte chemotactic protein 4, MCP-4, NCC-1, Small-inducible cytokine A13, MCP4, NCC1, SCYA13
Accession No.	Uniprot ID:Q99616
Uniprot	Q99616
GeneID	6357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: CCL13

Note: This product is for in vitro research use only