

C-C chemokine receptor type 5

Catalog No: #AP81457



Package Size: #AP81457-1 50ug #AP81457-2 100ug #AP81457-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-C CKR-5,MIP-1 alpha receptor,CD195,Cmkbr5
Accession No.	Uniprot ID:P51682
Uniprot	P51682
GeneID	12774;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccr5

Note: This product is for in vitro research use only