

B-lymphocyte antigen CD20

Catalog No: #AP81458



Package Size: #AP81458-1 50ug #AP81458-2 100ug #AP81458-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-lymphocyte antigen CD20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell differentiation antigen Ly-44,Lymphocyte antigen 44,Membrane-spanning 4-domains subfamily A member 1,CD20,Cd20,Ly-44,Ms4a2
Accession No.	Uniprot ID:P19437
Uniprot	P19437
GeneID	12482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ms4a1

Note: This product is for in vitro research use only