

Calcineurin subunit B type 1

Catalog No: #AP81461



Package Size: #AP81461-1 50ug #AP81461-2 100ug #AP81461-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcineurin subunit B type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Protein phosphatase 2B regulatory subunit 1,Protein phosphatase 3 regulatory subunit B alpha isoform 1,Cna2,Cnb
Accession No.	Uniprot ID:P63100
Uniprot	P63100
GeneID	29748;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppp3r1

Note: This product is for in vitro research use only