

# Cholecystokinin

Catalog No: #AP81465



Package Size: #AP81465-1 50ug #AP81465-2 100ug #AP81465-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Cholecystokinin                                        |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Rt                                                     |
| Other Names        | CCK                                                    |
| Accession No.      | Uniprot ID:P01355                                      |
| Uniprot            | P01355                                                 |
| GeneID             | 25298;                                                 |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cck

Note: This product is for in vitro research use only