

2-hydroxyacylsphingosine 1-beta-galactosyltransferase

Catalog No: #AP81466

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81466-1 50ug #AP81466-2 100ug #AP81466-3 1mg

Description

Product Name	2-hydroxyacylsphingosine 1-beta-galactosyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ceramide UDP-galactosyltransferase, Cerebroside synthase, UDP-galactose-ceramide galactosyltransferase, Cgt, Ugt4, Ugt8a
Accession No.	Uniprot ID:Q64676
Uniprot	Q64676
GeneID	22239;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ugt8

Note: This product is for in vitro research use only