

Caspase-4

Catalog No: #AP81468



Package Size: #AP81468-1 50ug #AP81468-2 100ug #AP81468-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CASP-4,Caspase-11,CASP-11,Protease ICH-3,Casp11,Caspl,lch3
Accession No.	Uniprot ID:P70343
Uniprot	P70343
GeneID	12363;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Casp4

Note: This product is for in vitro research use only