

Collagen type IV alpha-3-binding protein

Catalog No: #AP81473



Package Size: #AP81473-1 50ug #AP81473-2 100ug #AP81473-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen type IV alpha-3-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ceramide transfer protein,hCERT,Goodpasture antigen-binding protein,GPBP,START domain-containing protein 11,StARD11,StAR-related lipid transfer protein 11,CERT,STARD11
Accession No.	Uniprot ID:Q9Y5P4
Uniprot	Q9Y5P4
GeneID	10087;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COL4A3BP

Note: This product is for in vitro research use only