

Voltage-dependent P/Q-type calcium channel subunit alpha-1A

Catalog No: #AP81474

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81474-1 50ug #AP81474-2 100ug #AP81474-3 1mg

Description

Product Name	Voltage-dependent P/Q-type calcium channel subunit alpha-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain calcium channel I, BI, Calcium channel, L type, alpha-1 polypeptide isoform 4, Voltage-gated calcium channel subunit alpha Cav2.1, Caca1a, Cach4, Cacn3, Cacn1a4, Ccha1a
Accession No.	Uniprot ID:P97445
Uniprot	P97445
GeneID	12286;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Caca1a

Note: This product is for in vitro research use only