

Dipeptidyl peptidase 1

Catalog No: #AP81480



Package Size: #AP81480-1 50ug #AP81480-2 100ug #AP81480-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cathepsin C,Cathepsin J,Dipeptidyl peptidase I,DPP-I,Dipeptidyl transferase
Accession No.	Uniprot ID:P97821
Uniprot	P97821
GeneID	13032;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ctsc

Note: This product is for in vitro research use only