

Cathepsin K

Catalog No: #AP81482



Package Size: #AP81482-1 50ug #AP81482-2 100ug #AP81482-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cathepsin K
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Accession No.	Uniprot ID:Q3ZKN1
Uniprot	Q3ZKN1
GeneID	608843;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTSK

Note: This product is for in vitro research use only