

Bile salt-activated lipase

Catalog No: #AP81484



Package Size: #AP81484-1 50ug #AP81484-2 100ug #AP81484-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Bile salt-activated lipase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | BAL, Bile salt-stimulated lipase, BSSL, Carboxyl ester lipase, Cholesterol esterase, Pancreatic lysophospholipase, Sterol esterase, Lip1 |
| Accession No. | Uniprot ID: Q64285 |
| Uniprot | Q64285 |
| GeneID | 12613; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cel

Note: This product is for in vitro research use only