

C-C motif chemokine 5

Catalog No: #AP81488



Package Size: #AP81488-1 50ug #AP81488-2 100ug #AP81488-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Small-inducible cytokine A5,T-cell-specific protein RANTES,SCYA5
Accession No.	Uniprot ID:Q8HYS0
Uniprot	Q8HYS0
GeneID	403522;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL5

Note: This product is for in vitro research use only