

## C-C chemokine receptor type 5

Catalog No: #AP81489

Package Size: #AP81489-1 50ug #AP81489-2 100ug #AP81489-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C chemokine receptor type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-C CKR-5,MIP-1 alpha receptor,CD195,Cmkbr5
Accession No.	Uniprot ID:O08556
Uniprot	O08556
GeneID	117029;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ccr5

Note: This product is for in vitro research use only