

Caprin-1

Catalog No: #AP81491



Package Size: #AP81491-1 50ug #AP81491-2 100ug #AP81491-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caprin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytoplasmic activation- and proliferation-associated protein 1,GPI-anchored membrane protein 1,GPI-anchored protein p137,GPI-p137,Membrane component chromosome 11 surface marker 1,RNA granule protein 105,Gpiap,Gpiap1,Gpip137,M11s1,Rng105,G5E5
Accession No.	Uniprot ID:Q60865
Uniprot	Q60865
GeneID	53872;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Caprin1

Note: This product is for in vitro research use only